

## Anti-Mouse CD122/IL-2R $\beta$ Antibody (TM-Beta 1)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	The Anti-Mouse CD122/IL-2R $\beta$ Antibody (TM-Beta 1) is a murine-derived antibody of the IgG2bkappa class, specifically targeting CD122/IL-2R $\beta$ . Its isotype control is the Mouse IgG2bkappa Isotype Control.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481