

## HSD17B13 degrader 2

## Chemical Properties

CAS No. : 3035000-69-9


Formula: C<sub>21</sub>H<sub>19</sub>ClF<sub>3</sub>N<sub>3</sub>O<sub>2</sub>

Molecular Weight: 437.84

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	HSD17B13 degrader 2, a PROTAC target protein ligand (Ligands for Target Protein for PROTACs), is designed for synthesis applications.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2839 mL	11.4197 mL	22.8394 mL
5 mM	0.4568 mL	2.2839 mL	4.5679 mL
10 mM	0.2284 mL	1.142 mL	2.2839 mL
50 mM	0.0457 mL	0.2284 mL	0.4568 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481