

Anti-TIM-1/KIM-1/CD365 Antibody (2.70.2)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-TIM-1/KIM-1/CD365 Antibody (2.70.2) is a type of antibody targeting TIM/KIM-1/CD365.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481