


## (S,R,S)-AHPC-phenylacetic acid

### Chemical Properties

CAS No. :	2010985-19-8
Formula:	C30H36N4O4S
Molecular Weight:	548.7
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	(S,R,S)-AHPC-phenylacetic acid, an E3 Ligase Ligand-Linker Conjugate, facilitates the synthesis of PROTAC SMARCA2/4-degrader-15.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8225 mL	9.1124 mL	18.2249 mL
5 mM	0.3645 mL	1.8225 mL	3.645 mL
10 mM	0.1822 mL	0.9112 mL	1.8225 mL
50 mM	0.0364 mL	0.1822 mL	0.3645 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481