

## Anti-Mouse PD-L1/B7-H1 Antibody (10F.2H11)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	The Anti-Mouse PD-L1/B7-H1 Antibody (10F.2H11) is a rat IgG2b kappa mouse antibody specifically targeting PD-L1/B7-H1. The recommended isotype control for this antibody is Rat IgG2b kappa, Isotype Control.
Targets(IC50)	PD-1/PD-L1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481