

## E3 Ligase Ligand-linker Conjugate 112

## Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>34</sub>N<sub>4</sub>O<sub>3</sub>

Molecular Weight: 438.56

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	E3 Ligase Ligand-linker Conjugate 112, a type of E3 ubiquitinase ligand-linker conjugate, serves as a precursor for the synthesis of SMARCA2 degrader-10.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2802 mL	11.4009 mL	22.8019 mL
5 mM	0.456 mL	2.2802 mL	4.5604 mL
10 mM	0.228 mL	1.1401 mL	2.2802 mL
50 mM	0.0456 mL	0.228 mL	0.456 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481