

Anti-Mouse NKG2D/CD314 Antibody (HMG2D)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-Mouse NKG2D/CD314 Antibody (HMG2D) is an antibody of mouse origin, specifically an IgG derived from American hamster, targeting NKG2D/CD314. The isotype control for Anti-Mouse NKG2D/CD314 Antibody (HMG2D) is Armenian hamster IgG, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481