


## Arabinosylhypoxanthine

## Chemical Properties

CAS No. :	7013-16-3
Formula:	C10H12N4O5
Molecular Weight:	268.23
Storage:	Keep away from moisture, Store at low temperature Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Arabinosylhypoxanthine inhibits Ns5b polymerase in hepatitis C virus and can be used in biochemical experiments and organic synthesis.
Targets(IC50)	Nucleoside Antimetabolite/Analog, HCV Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7281 mL	18.6407 mL	37.2814 mL
5 mM	0.7456 mL	3.7281 mL	7.4563 mL
10 mM	0.3728 mL	1.8641 mL	3.7281 mL
50 mM	0.0746 mL	0.3728 mL	0.7456 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Robak T, Robak P. Purine nucleoside analogs in the treatment of rarer chronic lymphoid leukemias. *Curr Pharm Des.* 2012;18(23):3373-88.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481