


13,14-dihydro-15-keto Prostaglandin E1

Chemical Properties

CAS No. :	5094-14-4
Formula:	C ₂₀ H ₃₄ O ₅
Molecular Weight:	354.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	13,14-dihydro-15-keto Prostaglandin E1 inhibits ADP-induced platelet aggregation in human isolated platelet-rich plasma (IC ₅₀ =14.8 µg/mL) and is a PGE1 metabolite.
Targets(IC ₅₀)	Others, Drug Metabolite

Solubility Information

Solubility	PBS (pH 7.2): >1.6 mg/mL (from PGE1), Sonication is recommended. DMSO: >50 mg/mL (from PGE1), Sonication is recommended. DMF: >100 mg/mL (from PGE1), Sonication is recommended. Ethanol: >50 mg/mL (from PGE1), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.821 mL	14.1052 mL	28.2103 mL
5 mM	0.5642 mL	2.821 mL	5.6421 mL
10 mM	0.2821 mL	1.4105 mL	2.821 mL
50 mM	0.0564 mL	0.2821 mL	0.5642 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481