


Polygalasaponin E

Chemical Properties

CAS No. : 882664-72-4
 Formula: C₄₇H₇₄O₁₉
 Molecular Weight: 943.09
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Polygalasaponin E is a natural product for research related to life sciences. The catalog number is TN6531 and the CAS number is 882664-72-4.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0603 mL	5.3017 mL	10.6034 mL
5 mM	0.2121 mL	1.0603 mL	2.1207 mL
10 mM	0.106 mL	0.5302 mL	1.0603 mL
50 mM	0.0212 mL	0.106 mL	0.2121 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481