

BIMAX1 acetate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

BIMAX1 acetate is a tool peptide utilized for investigating the interactions between RBBP4 and nuclear import protein (Importin) β 1, as well as their impact on the nuclear import pathway. Additionally, BIMAX1 acetate serves in studying the role of RBBP4 in regulating nuclear import efficiency and cellular senescence.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481