

PRMT3-IN-4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PRMT3-IN-4 (intermediate 15) is an inhibitor of Protein arginine methyltransferase 3 (PRMT3) and serves as the active control for SGC707. It can be utilized in the synthesis of PROTACs targeting PRMT3 and is applicable in research related to leukemia.
Targets(IC50)	Histone Methyltransferase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481