


## PROTAC SMARCA2/4-degrader-32

## Chemical Properties

CAS No. :	2568273-79-8
Formula:	C54H69N13O5S
Molecular Weight:	1012.28
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	PROTAC SMARCA2/4-degrader-32 (compound I-446) is a PROTAC degrader targeting both SMARCA2 and SMARCA4 proteins. In A549 cells, it shows DC50s values for the degradation of SMARCA2/4 proteins 100 nM and achieves a maximum degradation rate (Dmax%) of >90 after 24 hours of treatment.
Targets(IC50)	PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9879 mL	4.9393 mL	9.8787 mL
5 mM	0.1976 mL	0.9879 mL	1.9757 mL
10 mM	0.0988 mL	0.4939 mL	0.9879 mL
50 mM	0.0198 mL	0.0988 mL	0.1976 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481