

BI1543673

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BI1543673 is an inhibitor of interleukin-1 receptor-associated kinase 4 (IRAK4). It can diminish inflammatory responses in human lung tissue stimulated by TLR4 and TLR7/8. Additionally, BI1543673 reduces inflammatory signaling in a mouse lung inflammation model induced by LPS.
Targets(IC50)	IRAK,TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481