

SJ-106C

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	SJ-106C, a SIRT inhibitor, exhibits IC50 values of 0.59, 0.12, and 0.49 $\mu$ M against SIRT1/2/3, respectively. This compound specifically targets mitochondria and inhibits the growth of DLBCL tumors.
Targets(IC50)	Mitochondrial Metabolism,Sirtuin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481