

PD-1/PD-L1-IN-51

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PD-1/PD-L1-IN-51 (Compound III-4) is an inhibitor of PD-1/PD-L1 (IC50: hPD-L1 at 2.9 nM). It binds directly to PD-L1, blocking the interaction between PD-1 and PD-L1, and enhancing the release of IFN- γ . Additionally, PD-1/PD-L1-IN-51 exhibits antitumor activity.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481