

## PROTAC SMARCA2/4-degrader-22

### Chemical Properties

CAS No. : 2568277-48-3


Formula: C<sub>51</sub>H<sub>61</sub>N<sub>11</sub>O<sub>6</sub>S

Molecular Weight: 956.17

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	SMARCA2 degrader-14 (compound I-408) is a PROTAC degrader targeting both SMARCA2 and SMARCA4. In A549 cells, it achieves degradation of SMARCA2/4 proteins with DC50s less than 100 nM and a maximal degradation rate (Dmax%) exceeding 90% after a 24-hour treatment.
Targets(IC50)	PROTACs

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0458 mL	5.2292 mL	10.4584 mL
5 mM	0.2092 mL	1.0458 mL	2.0917 mL
10 mM	0.1046 mL	0.5229 mL	1.0458 mL
50 mM	0.0209 mL	0.1046 mL	0.2092 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481