


## PROTAC SMARCA2/4-degrader-33

### Chemical Properties

CAS No. :	2568276-36-6
Formula:	C44H51N11O4
Molecular Weight:	797.95
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	PROTAC SMARCA2/4-degrader-33 (compound I-277) functions as a degrader targeting the proteins SMARCA2 and SMARCA4, exhibiting DC50s values of less than 100 nM for the degradation of SMARCA2/4 in A549 cells. After a treatment period of 24 hours, the maximum degradation rates (Dmax%) consistently exceed 90%.
Targets(IC50)	PROTACs

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2532 mL	6.2661 mL	12.5321 mL
5 mM	0.2506 mL	1.2532 mL	2.5064 mL
10 mM	0.1253 mL	0.6266 mL	1.2532 mL
50 mM	0.0251 mL	0.1253 mL	0.2506 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481