

Antifeedant agent 1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Antifeedant agent 1 functions as a carboxylesterase inhibitor that controls insect pests by inhibiting their feeding, with an EC50 of 0.038 mg/mL. It is applicable in the research of pest control.
Targets(IC50)	Lipid

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481