


## PROTAC SMARCA2 degrader-10

## Chemical Properties

CAS No. :	2568279-34-3
Formula:	C55H60N10O5S
Molecular Weight:	973.19
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	PROTAC SMARCA2 degrader-10 (compound I-507) is a targeted degrader of SMARCA2, demonstrating a DC50 of less than 100 nM in A549 cells, with a maximum degradation rate (Dmax%) exceeding 90% after 24 hours of treatment.
Targets(IC50)	PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0275 mL	5.1377 mL	10.2755 mL
5 mM	0.2055 mL	1.0275 mL	2.0551 mL
10 mM	0.1028 mL	0.5138 mL	1.0275 mL
50 mM	0.0206 mL	0.1028 mL	0.2055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481