

Dansyl-Glu-Gly-Arg-Chloromethylketone TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dansyl-Glu-Gly-Arg-Chloromethylketone TFA is a protease inhibitor that is effective against serine/threonine proteases. It also inhibits activated porcine factor IX.
Targets(IC50)	Serine Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481