

KL-465

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	KL-465, a MAGE-A3 PROTAC degrader, effectively degrades MAGE-A3 in HeLa cells with a DC50 of 2 μ M. It also inhibits the proliferation of MAGE-A3-positive cancer cells, specifically in HCT116, A375, and HeLa cell lines.
Targets(IC50)	PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481