

TRPV1/CB2 agonist 1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TRPV1/CB2 agonist 1 (compound 41) is a dual agonist targeting hTRPV1 and CB2 receptors, with an EC50 value of 26.8 μ M for hTRPV1. It is applicable in research related to the nervous system.
Targets(IC50)	Cannabinoid Receptor,TRP/TRPV Channel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481