

## Antifungal agent 110

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Antifungalagent 110 exhibits antifungal activity against Botrytis cinerea and Rhizoctonia solani. It demonstrates significant effectiveness in inhibiting R. solani, with an EC50 value of 0.27 mg/L. Antifungalagent 110 shows promise for use in plant disease research.
Targets(IC50)	Antifungal

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481