


F1386-0303

## Chemical Properties

CAS No. : 287177-12-2  
Formula: C<sub>18</sub>H<sub>13</sub>N<sub>3</sub>O  
Molecular Weight: 287.32  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	F1386-0303 is a novel potent and selective MAP4K4 inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4804 mL	17.4022 mL	34.8044 mL
5 mM	0.6961 mL	3.4804 mL	6.9609 mL
10 mM	0.348 mL	1.7402 mL	3.4804 mL
50 mM	0.0696 mL	0.348 mL	0.6961 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481