

## Eyeseryl

## Chemical Properties

CAS No. : 820959-17-9


Formula: C<sub>20</sub>H<sub>28</sub>N<sub>8</sub>O<sub>7</sub>

Molecular Weight: 492.49

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Eyeseryl, also known as Acetyl Tetrapeptide-5, humectant or hygroscopic moisturizer which helps reduce eye puffiness, improve skin elasticity as well as overall smoothness.
Targets(IC <sub>50</sub> )	Angiotensin-converting Enzyme (ACE)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0305 mL	10.1525 mL	20.305 mL
5 mM	0.4061 mL	2.0305 mL	4.061 mL
10 mM	0.203 mL	1.0152 mL	2.0305 mL
50 mM	0.0406 mL	0.203 mL	0.4061 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481