

## Eurycomalin A

## Chemical Properties


CAS No. :

Formula: C<sub>24</sub>H<sub>22</sub>O<sub>6</sub>

Molecular Weight: 406.434

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Eurycomalin A is a useful organic compound for research related to life sciences and the catalog number is T124210.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4604 mL	12.3022 mL	24.6045 mL
5 mM	0.4921 mL	2.4604 mL	4.9209 mL
10 mM	0.246 mL	1.2302 mL	2.4604 mL
50 mM	0.0492 mL	0.246 mL	0.4921 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481