

Cemtirestat

Chemical Properties

| | |
|-------------------|--|
| CAS No. : | 309283-89-4 |
| Formula: | C11H8N4O2S |
| Molecular Weight: | 260.272 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small> |

Biological Description

| | |
|---------------|--|
| Description | Cemtirestat is an aldose reductase inhibitor |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 3.8422 mL | 19.2108 mL | 38.4216 mL |
| 5 mM | 0.7684 mL | 3.8422 mL | 7.6843 mL |
| 10 mM | 0.3842 mL | 1.9211 mL | 3.8422 mL |
| 50 mM | 0.0768 mL | 0.3842 mL | 0.7684 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481