

Phenolic acid

Chemical Properties


CAS No. : 29656-58-4

Formula: C7H8D2O3

Molecular Weight: 144.16

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenolic acid is an endogenous metabolite.
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9367 mL	34.6837 mL	69.3674 mL
5 mM	1.3873 mL	6.9367 mL	13.8735 mL
10 mM	0.6937 mL	3.4684 mL	6.9367 mL
50 mM	0.1387 mL	0.6937 mL	1.3873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481