


## (S)-GNA-U-phosphoramidite

## Chemical Properties

CAS No. :	494784-15-5
Formula:	C37H45N4O7P
Molecular Weight:	688.75
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-GNA-U-phosphoramidite is a Nucleoside Phosphoramidite.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4519 mL	7.2595 mL	14.5191 mL
5 mM	0.2904 mL	1.4519 mL	2.9038 mL
10 mM	0.1452 mL	0.726 mL	1.4519 mL
50 mM	0.029 mL	0.1452 mL	0.2904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481