

Isoorientin 2''-O-rhamnoside

Chemical Properties


CAS No. :

Formula: C54H58O30

Molecular Weight: 1187.03

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isoorientin 2''-O-rhamnoside is a useful organic compound for research related to life sciences and the catalog number is T124742.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8424 mL	4.2122 mL	8.4244 mL
5 mM	0.1685 mL	0.8424 mL	1.6849 mL
10 mM	0.0842 mL	0.4212 mL	0.8424 mL
50 mM	0.0168 mL	0.0842 mL	0.1685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481