


## Osimertinib

## Chemical Properties

CAS No. :	1421373-65-0
Formula:	C <sub>28</sub> H <sub>33</sub> N <sub>7</sub> O <sub>2</sub>
Molecular Weight:	499.61
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Osimertinib (AZD-9291) is an EGFR third-generation inhibitor that inhibits the T790M resistance mutation produced by second-generation EGFR inhibitors with irreversible and oral activity. Osimertinib has antitumor activity for the treatment of EGFR-mutated non-small-cell lung cancer.
Targets(IC50)	EGFR
In vitro	In the EGFRm+ (PC9) and EGFRm+/T790M (H1975) tumor models, oral administration of AZD9291 (5 mg/kg) effectively inhibits AKT and ERK signaling pathways as well as EGFR phosphorylation in tumors, leading to tumor regression.
In vivo	In mutant EGFR cell lines, AZD9291 effectively inhibits cell proliferation.
Kinase Assay	ACY-1215 is dissolved and subsequently diluted in assay buffer [50 mM HEPES, pH 7.4, 100 mM KCl, 0.001% Tween-20, 0.05% BSA, and 20 μM tris(2-carboxyethyl)phosphine] to 6-fold the final concentration. HDAC enzymes are diluted to 1.5-fold of the final concentration in assay buffer and pre-incubated with ACY-1215 for 10 minutes before the addition of the substrate. The amount of FTS (HDAC1, HDAC2, HDAC3, and HDAC6) or MAZ-1675 (HDAC4, HDAC5, HDAC7, HDAC8, and HDAC9) used for each enzyme is equal to the Michaelis constant (K <sub>m</sub> ), as determined by a titration curve. FTS or MAZ-1675 is diluted in assay buffer to 6-fold the final concentration with 0.3 μM sequencing grade trypsin. The substrate/trypsin mix is added to the enzyme/compound mix and the plate is shaken for 60 seconds and then placed into a SpectraMax M5 microtiter plate reader. The enzymatic reaction is monitored for release of 7-amino-4-methoxy-coumarin over 30 minutes, after deacetylation of the lysine side chain in the peptide substrate, and the linear rate of the reaction is calculated[1].
Cell Research	AZD-9291 is dissolved in DMSO and stored, and then diluted with appropriate medium before use[1]. PC-9 cells are seeded into T75 flasks (5×10 <sup>5</sup> cells/flask) in RPMI growth media and incubated at 37°C, 5% CO <sub>2</sub> . The following day the media is replaced with media supplemented with a concentration of EGFR inhibitor equal to the EC <sub>50</sub> concentration predetermined in PC-9 cells. Media changes are carried out every 2-3 days and resistant clones allowed to grow to 80% confluency prior to the cells being trypsinised and reseeded at the original seeding density in media containing twice the concentration of EGFR inhibitor. Dose escalations are continued until a final concentration of 1.5 μM Gefitinib, 1.5 μM Afatinib, 1.5 μM WZ4002 or 160 nM AZD-9291

## A DRUG SCREENING EXPERT

Cell Research	are achieved[1].
---------------	------------------

### Solubility Information

Solubility	Ethanol: 22 mg/mL (44.03 mM),Sonication is recommended. DMSO: 120 mg/mL (240.19 mM),Sonication is recommended. H2O: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 9.2 mg/mL (18.41 mM),Suspension. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0016 mL	10.0078 mL	20.0156 mL
5 mM	0.4003 mL	2.0016 mL	4.0031 mL
10 mM	0.2002 mL	1.0008 mL	2.0016 mL
50 mM	0.040 mL	0.2002 mL	0.4003 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

### Reference

- Hirano T, et al. Pharmacological and Structural Characterizations of Naquotinib, a Novel Third-Generation EGFR Tyrosine Kinase Inhibitor, in EGFR-Mutated Non-Small Cell Lung Cancer. *Mol Cancer Ther.* 2018 Apr;17(4):740-750.
- Liu Y, Lai M, Li S, et al. LS-106, a novel EGFR inhibitor targeting C797S, exhibits anti-tumor activities both in vitro and in vivo. *Cancer Science.* 2021
- Zhang C, Sun Y X, Yi D C, et al. Neoadjuvant sintilimab plus chemotherapy in EGFR-mutant NSCLC: Phase 2 trial interim results (NEOTIDE/CTONG2104). *Cell Reports Medicine.* 2024
- Wang J, Wang Y, Zhang S, et al. Inhibitory effect of 1, 4, 5, 6-tetrahydroxy-7, 8-diprenylxanthone against NSCLC with L858R/T790M/C797S mutant EGFR. *Scientific Reports.* 2024, 14(1): 26549.
- Chen Y, Jin Y, Wang W, et al. Preparation and antitumor activity of selenium nanocomposite stabilized by AAGL from *Agerocybe aegerita*. *International Journal of Biological Macromolecules.* 2024: 137002.
- Peng T, Ma X, Hua W, et al. Individualized patient tumor organoids faithfully preserve human brain tumor ecosystems and predict patient response to therapy. *Cell Stem Cell.* 2025
- Lv S, Pan Q, Lu W, et al. Tenovin 3 induces apoptosis and ferroptosis in EGFR 19del non small cell lung cancer cells. *Scientific Reports.* 2024, 14(1): 7654.
- Cross DA, et al. AZD9291, an irreversible EGFR TKI, overcomes T790M-mediated resistance to EGFR inhibitors in lung cancer. *Cancer Discov.* 2014 Sep;4(9):1046-61.
- Yu J, Zhang L, Peng J, et al. Dictamnine, a novel c-Met inhibitor, suppresses the proliferation of lung cancer cells by downregulating the PI3K/AKT/mTOR and MAPK signaling pathways. *Biochemical pharmacology.* 2022, 195: 114864.
- Wang N, et al. Osimertinib (AZD9291) increases radio-sensitivity in EGFR T790M non-small cell lung cancer. *Oncol Rep.* 2019 Jan;41(1):77-86.
- Finlay MR, et al. Discovery of a potent and selective EGFR inhibitor (AZD9291) of both sensitizing and T790M resistance mutations that spares the wild type form of the receptor. *J Med Chem.* 2014 Oct 23;57(20):8249-67.
- Shang J, Ning S, Chen Y, et al. MDL-800, an allosteric activator of SIRT6, suppresses proliferation and enhances EGFR-TKIs therapy in non-small cell lung cancer. *Acta Pharmacologica Sinica.* 2021, 42(1): 120-131
- Jin Y, Zhang Y, Huang A, et al. Overexpression of SERPINA3 suppresses tumor progression by modulating SPOP/NF- $\kappa$ B in lung cancer. *International Journal of Oncology.* 2023, 63(2): 1-16.
- Nanjo S, et al. High efficacy of third generation EGFR inhibitor AZD9291 in a leptomeningeal carcinomatosis model with EGFR-mutant lung cancer cells. *Oncotarget.* 2016 Jan 26;7(4):3847-56.
- Liang J, Bi G, Sui Q, et al. Transcription factor ZNF263 enhances EGFR-targeted therapeutic response and reduces residual disease in lung adenocarcinoma. *Cell Reports.* 2024, 43(2).
- Tan J, Zhu L, Shi J, et al. Evaluation of drug resistance for EGFR-TKIs in Lung Cancer via Multicellular lung-on-a-chip. *European Journal of Pharmaceutical Sciences.* 2024: 106805.
- Zhang Y, Chen Y, Wang B, et al. A novel selenium nanocomposite modified by AANL inhibits tumor growth by upregulating CLK2 in lung cancer. *Bioorganic Chemistry.* 2024: 107459.
- Wei X, Zhang G, Liu Q, et al. Almonertinib and alflutinib show novel inhibition on rare EGFR S768I mutant cells. *Clinical and Translational Oncology.* 2024: 1-16.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel: 781-999-4286    E\_mail: info@targetmol.com    Address: 34 Washington Street, Wellesley Hills, MA 02481