


PYR-0503

Chemical Properties

CAS No. : 16200-50-3
 Formula: C₁₂H₁₉NO₂
 Molecular Weight: 209.28
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	PYR-0503 is an intermediate used for chemical synthesis of a number of biologically important molecules, including bile pigments, photosensitizers, porphyrins, anticancer drugs.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, H ₂ O: Insoluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7783 mL	23.8914 mL	47.7829 mL
5 mM	0.9557 mL	4.7783 mL	9.5566 mL
10 mM	0.4778 mL	2.3891 mL	4.7783 mL
50 mM	0.0956 mL	0.4778 mL	0.9557 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481