

Antifungal agent 116

Chemical Properties

CAS No. :

Formula: C₁₂H₁₀FN₃

Molecular Weight: 215.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Antifungalagent 116 (compound 1a) exhibits broad-spectrum antifungal (fungal) activity.
Targets(IC50)	Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6462 mL	23.231 mL	46.4619 mL
5 mM	0.9292 mL	4.6462 mL	9.2924 mL
10 mM	0.4646 mL	2.3231 mL	4.6462 mL
50 mM	0.0929 mL	0.4646 mL	0.9292 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481