

ERAP1-IN-2

Chemical Properties

CAS No. :

Formula: C₁₁H₉NO₃S₂

Molecular Weight: 267.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ERAP1-IN-2 (compound 3f) serves as an inhibitor for ERAP1 (endoplasmic reticulum aminopeptidase 1). The IC ₅₀ value of ERAP1-IN-2 against hERAP1 is 1.72 μM.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7408 mL	18.7042 mL	37.4083 mL
5 mM	0.7482 mL	3.7408 mL	7.4817 mL
10 mM	0.3741 mL	1.8704 mL	3.7408 mL
50 mM	0.0748 mL	0.3741 mL	0.7482 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481