


EV206

Chemical Properties

CAS No. : 2247047-81-8
 Formula: C₂₁H₁₉N₃O
 Molecular Weight: 329.40
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	EV206 is an effective inhibitor of Hsp70, exhibiting antiproliferative properties. This compound induces cell apoptosis (apoptosis) and enhances the protein expression of cleaved PARP and caspase-3. Additionally, EV206 demonstrates antitumor activity.
Targets(IC50)	HSP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0358 mL	15.1791 mL	30.3582 mL
5 mM	0.6072 mL	3.0358 mL	6.0716 mL
10 mM	0.3036 mL	1.5179 mL	3.0358 mL
50 mM	0.0607 mL	0.3036 mL	0.6072 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481