

PROTAC SMARCA2 degrader-18

Chemical Properties

CAS No. : 2408393-47-3


Formula: C₄₀H₄₁N₉O₃

Molecular Weight: 695.81

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PROTAC SMARCA2 degrader-18 (Example144) is a PROTAC SMARCA2 degrader with potential for researching non-small cell lung cancer (NSCLC).
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4372 mL	7.1859 mL	14.3717 mL
5 mM	0.2874 mL	1.4372 mL	2.8743 mL
10 mM	0.1437 mL	0.7186 mL	1.4372 mL
50 mM	0.0287 mL	0.1437 mL	0.2874 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481