

PROTAC BRM/BRG1 degrader-4

Chemical Properties

CAS No. :

Formula: C44H52N10O6S

Molecular Weight: 849.01

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PROTAC BRM/BRG1 degrader-4 (example 047) is a potent degrader targeting both BRM and BRG1, demonstrating DC50 values of 0.21 nM and 3.4 nM respectively in SW1573 cells.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1778 mL	5.8892 mL	11.7784 mL
5 mM	0.2356 mL	1.1778 mL	2.3557 mL
10 mM	0.1178 mL	0.5889 mL	1.1778 mL
50 mM	0.0236 mL	0.1178 mL	0.2356 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481