


CB1/2 receptor-1

Chemical Properties

CAS No. : 1260669-31-5
 Formula: C33H48O2
 Molecular Weight: 476.73
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	CB1/2 receptor-1 (compound 5.3) serves as an agonist for CB1/2 receptors and is utilized in angiogenesis research.
Targets(IC50)	Cannabinoid Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0976 mL	10.4881 mL	20.9762 mL
5 mM	0.4195 mL	2.0976 mL	4.1952 mL
10 mM	0.2098 mL	1.0488 mL	2.0976 mL
50 mM	0.042 mL	0.2098 mL	0.4195 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481