

## Acaricide

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>21</sub>F<sub>2</sub>NOS<sub>2</sub>

Molecular Weight: 405.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Compound 7l, an acaricide, serves as an effective LSD1 inhibitor. It demonstrates substantial acaricidal activity against the eggs of <i>T. cinnabarinus</i> , with a lethal concentration 50 (LC50) of 0.0035 mg/L.
Targets(IC50)	Histone Demethylase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.466 mL	12.3298 mL	24.6597 mL
5 mM	0.4932 mL	2.466 mL	4.9319 mL
10 mM	0.2466 mL	1.233 mL	2.466 mL
50 mM	0.0493 mL	0.2466 mL	0.4932 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481