

YAP/TAZ inhibitor-4

Chemical Properties

CAS No. :

Formula: C₂₈H₂₈F₃N₃O₃

Molecular Weight: 511.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	YAP/TAZ Inhibitor-4 (Compound 45) acts as an inhibitor of YAP/TAZ, exhibiting inhibitory activity on TEAD transcriptional activation by preventing the interaction between YAP or TAZ and TEAD.
Targets(IC50)	YAP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9549 mL	9.7744 mL	19.5488 mL
5 mM	0.391 mL	1.9549 mL	3.9098 mL
10 mM	0.1955 mL	0.9774 mL	1.9549 mL
50 mM	0.0391 mL	0.1955 mL	0.391 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481