

## Pediocin PA 1 TFA

### Chemical Properties

CAS No. :

Formula: C196H293N61O60S5.xC2HF3O2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Pediocin PA 1 TFA is a broad-spectrum bacteriocin produced by lactic acid bacteria, inhibiting the activity of foodborne pathogens such as <i>Listeria monocytogenes</i> and other Gram-positive bacteria. It has applications as a biopreservative in food products.
Targets(IC50)	Antibacterial

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481