

Tyk2-IN-20

## Chemical Properties


CAS No. : 3058568-08-1

Formula: C<sub>24</sub>H<sub>25</sub>N<sub>7</sub>O<sub>2</sub>

Molecular Weight: 443.50

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tyk2-IN-20 (Example 289) is an effective inhibitor of Tyk2 with an IC <sub>50</sub> value below 5 nM. Additionally, it inhibits JAK1, JAK2, and JAK3 with IC <sub>50</sub> values under 100 nM. This compound is utilized for the research of inflammatory diseases.
Targets(IC <sub>50</sub> )	JAK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2548 mL	11.274 mL	22.5479 mL
5 mM	0.451 mL	2.2548 mL	4.5096 mL
10 mM	0.2255 mL	1.1274 mL	2.2548 mL
50 mM	0.0451 mL	0.2255 mL	0.451 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481