

Multi-target kinase inhibitor 3

Chemical Properties

CAS No. :

Formula: C₂₃H₂₂N₆O₂S

Molecular Weight: 446.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Multi-target kinase inhibitor 3 (compound 6i) is a potent multi-target anticancer agent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2395 mL	11.1977 mL	22.3954 mL
5 mM	0.4479 mL	2.2395 mL	4.4791 mL
10 mM	0.224 mL	1.1198 mL	2.2395 mL
50 mM	0.0448 mL	0.224 mL	0.4479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481