

sPLA2 inhibitor 2

Chemical Properties

CAS No. :

Formula: C₂₀H₁₈N₆O₄S₂

Molecular Weight: 470.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	sPLA2 inhibitor 2 (compound 6a) is an sPLA2 inhibitor with an IC ₅₀ value of 0.0475 μM, making it valuable for diabetes research.
Targets(IC ₅₀)	Phospholipase,Virus Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1253 mL	10.6265 mL	21.2531 mL
5 mM	0.4251 mL	2.1253 mL	4.2506 mL
10 mM	0.2125 mL	1.0627 mL	2.1253 mL
50 mM	0.0425 mL	0.2125 mL	0.4251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481