


## Androgen receptor antagonist 11

## Chemical Properties

CAS No. :	2719816-32-5
Formula:	C <sub>20</sub> H <sub>19</sub> F <sub>3</sub> N <sub>4</sub> O <sub>3</sub> S
Molecular Weight:	452.45
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Androgen receptor antagonist 11 (compound N29) is a selective, orally available antagonist.
Targets(IC50)	Others,Androgen Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2102 mL	11.0509 mL	22.1019 mL
5 mM	0.442 mL	2.2102 mL	4.4204 mL
10 mM	0.221 mL	1.1051 mL	2.2102 mL
50 mM	0.0442 mL	0.221 mL	0.442 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481