


(2R,3R)-2,3,4-Trihydroxybutanoic acid

Chemical Properties

CAS No. :	488-16-4
Formula:	C ₄ H ₈ O ₅
Molecular Weight:	136.1
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2R,3R)-2,3,4-Trihydroxybutanoic acid, also known as erythorbic acid, is a naturally occurring organic acid found in various bodily fluids including aqueous humor, urine, plasma, cerebrospinal fluid, and synovial fluid. It is prevalent in both healthy adults and children and is produced through the oxidation of GlcNAc.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.3475 mL	36.7377 mL	73.4754 mL
5 mM	1.4695 mL	7.3475 mL	14.6951 mL
10 mM	0.7348 mL	3.6738 mL	7.3475 mL
50 mM	0.147 mL	0.7348 mL	1.4695 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481