

PROTAC SMARCA2/4-degrader-35

Chemical Properties

CAS No. : 2933253-85-9


Formula: C56H69N13O6

Molecular Weight: 1020.23

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PROTAC SMARCA2/4-degrader-35 (Ex.43) is a degrader of SMARCA2/4 with a DC50 value of less than 2.5 nM.
Targets(IC50)	PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9802 mL	4.9009 mL	9.8017 mL
5 mM	0.196 mL	0.9802 mL	1.9603 mL
10 mM	0.098 mL	0.4901 mL	0.9802 mL
50 mM	0.0196 mL	0.098 mL	0.196 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481