


## PROTAC SMARCA2 degrader-23

## Chemical Properties

CAS No. :	2892523-67-8
Formula:	C47H57N9O6S
Molecular Weight:	876.08
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	PROTAC SMARCA2 degrader-23 (Example 1) is an effective and selective degrader of PROTAC SMARCA2, with a DC50 value of less than 100 nM. It has potential for use in cancer research.
Targets(IC50)	Epigenetic Reader Domain, PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1414 mL	5.7072 mL	11.4145 mL
5 mM	0.2283 mL	1.1414 mL	2.2829 mL
10 mM	0.1141 mL	0.5707 mL	1.1414 mL
50 mM	0.0228 mL	0.1141 mL	0.2283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481