


## Disperse orange 3

## Chemical Properties

CAS No. :	730-40-5
Formula:	C <sub>12</sub> H <sub>10</sub> N <sub>4</sub> O <sub>2</sub>
Molecular Weight:	242.23
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Disperse Orange 3 is a monoazo dye characterized by a dye content of 90%, along with dispersants and surfactants. It features two aromatic rings, one with an amino group and the other with a nitro group, displaying a distinct absorption peak at 415 nm.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1283 mL	20.6415 mL	41.2831 mL
5 mM	0.8257 mL	4.1283 mL	8.2566 mL
10 mM	0.4128 mL	2.0642 mL	4.1283 mL
50 mM	0.0826 mL	0.4128 mL	0.8257 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481